

# LIBS Sci-Trace

## Configurable spectroscopic system

Possibility of selecting desired combination of the interaction chamber, laser, detection system, specialized module and optomechanical accessories. Mutual compatibility guaranteed.

## Plug-in concept

Easy system expansion by the user, wide spectrum of modules: components of the interaction chamber, lasers, spectrometers, detectors, vacuum components etc. The user simply joins the new module into the setup and activates the corresponding software plug-in.

Sci-Trace is a laboratory LIBS research setup consisting of an instrumentation cabinet and the LIBS Interaction Chamber mounted on an optical breadboard.

Cabinet provides space for multiple shelfs with LIBS instrumentation (laser head, spectrometers, power meters, calibration lamps, etc.) and rack-compatible components (control electronics, PC, laser PSU etc.).

## Designed by scientists for scientists

Designed to be opened and ready for various researcher's extensions and experiments. Allows to fully concentrate on the LIBS method and its results rather than troubleshooting and system building.

### Capable and intuitive software

Integrated software for control of all the system elements, spectra capturing and spectra processing with still growing chemometric capabilities.





## LIBS Sci-Trace

#### Sample area

Manipulator		Motorized manipulator Movement range 60×80×45 mm, 2 µm resolution, High-Vacuum ready, including series of sample holders
Interaction chamber		Vacuum LIBS Interaction Chamber Airtight rigid body, 11 input ports aiming to a common center + 4 lateral input ports
[a	-	Open type - Cage chamber 6 mounting ports aiming to a common center
[8	alt]	Manipulator stand
Top optical breadboard		Anodized Al-alloy board, M6 threaded holes, dimensions: 1304 × 829 × 8 mm Feedthroughs for cable management and laser beam delivery, USB connector panel, possible to mount corner rails for laser filter plates
[a	alt]	Imperial threaded holes, magnetic steel board and different board dimensions upon request

#### **LIBS instruments**

[alt] alternative configuration [opt] optional feature

Pulsed laser	Lamp-pumped (LPSS) Nd:YAG 532nm, 200 mJ, 8 ns, 20 Hz, compact design
[a	lt] Diode-pumped Nd:YAG (DPSS)
[0]	ot] Double pulsed feature

[opt] Other Nd:YAG wavelengths available (1064nm, 532 nm, 355nm, 266 nm) [opt] Up to 800 mJ at 1064 (FPSS, single-pulsed)

SpectrometerEchelle, 190-1100 nmFocal length 120 nm, f/4, resolving power up to 5000 λ/FWHM

[alt] Czerny-Turner, multiple gratings on turret, USB control, multiple output

Detector EMCCD, 180-1100 nm

1004 x 1002 px, 20 Hz, min 10 µs exposuretime
[alt] iCCD detector, 1024 x 1024 px, 180-850 nm, USB
[alt] Deep-UV CCD detector (Chamber mounted)

Digital Delay Generator 4 output, 5 ns time resolution


[alt] 8 output, 5 ns time resolution

Accessories Calibration lamp - continuous spectrum: Deuterium-Halogen

Calibration lamp - line spectrum: Mercury-Argon

Guiding laser, DPSS 532 nm, 4 mW  $\,$ 

Laser safety glasses, 35% visible light transmission, OD 7+ (190-534 nm), OD 6+ (925-1070 nm)



### **LIBS Sci-Trace**

#### **Extension modules**

Pressure regulating system (PRS) Setting the pressure in the chamber in the range 1-1300 mbar (a), Ar, He, CO<sub>2</sub>, etc

Gas jet and exhaust module Cleaning the sample, creating local atmosphere of inert gas, pulse mode

Deep-UV CCD detector module Chamber port-mounted for detecting the emission lines in the region bellow 200 nm

Motorized defocusing module Chamber internal module for defocusing the laser beam by moving the focusing lens

Plasma imaging module Chamber internal module for defocusing the laser beam by moving the focusing le

Chamber internal module for defocusing the laser beam by moving the focusing le

Chamber port-mounted module with triggered CMOS camera (global shutter)

Magnetic field module Chamber internal module for confining the plasma in the magnetic field

**Liquid LIBS module** Chamber internal module for analyzing the liquids

Laser power meter module Realtime recording the laser energy value

Laser attenuator module Controling the laser energy while keeping the laser to operate at its most stable output power

#### Software capabilities

Manipulator movement

Sample view, laser autofocus

Chemical mapping, depth profiling

Control of connected LIBS instruments (lasers, detectors, DDGs, etc.)


Spectra capturing and manipulation

Identification of emission lines and chemical elements

Database of elements (spectra captured by LIBS)

Calculations of plasma parameters (temperature, electron density)

Creation of calibration curves



#### Instrumentation cabinet

Internal shelf system for LIBS instruments 2 anodized Al-alloy shelfs, M6 threaded holes

for mounting the laser head / spektroscope / optomechanics

[alt] Other number of shelfs upon request

Internal rack system for control electronics 19-inches rack, height 16U

Installed control electronics, control PC, laser PSU, DDG, PSR

Safety elements Interlock system on chamber door and cabinet door

Laser beam hidden in tubes

I/O panel 2x HDMI (dual monitor support), LAN, GAS inlet, Gas outlet, Vacuum pump output, Mains

Control panel Emergency STOP, key ON/OFF, electronics ON/OFF, PC ON/OFF, USB

**Housing and construction**Al profile frame covered by steel plates

Cooling fans, noise dumpening materials 4 doors: 2 for rack systam, 2 for shelf system

4 wheels with rectractable stands

**Dimensions & Weight** 1314 × 851 × 1471 mm, 330 kg (in the default configuration)

Power requirements ~230 V, 50 Hz, 16 A

