Distributed by:

soliton



PICAS Photonically Integratec Cold Atom Source

Related Products

The PICAS is used in conjunction with ColdQuanta Cold Atom Source Cells:

Fixed Source: CASC-2016-XX Exchangeable: CSJXX

Product

Alignment-free, 2D+ MOT source

Robust, fixed optics and opto-mechanics

Compact unit with flexible orientation

Fiber-coupled input

Integrated source cell

Available with rubidium or cesium

Flexible orientation



Product Description

ColdQuanta's PICAS, Photonically Integrated Cold Atom Source, is a compact, alignment-free package that produces a high-flux beam of laser-cooled atoms, using user-supplied, fiber-coupled light. When combined with a ColdQuanta source cell, the PICAS is easily integrated into the end user's vacuum system through a standard CF interface, and the opto-mechanics unit can be removed during bake-out.

The PICAS employs our patent-pending, ultra-compact beam expansion system. This unit produces a rectangular beam of circularly polarized light with an approximately uniform intensity profile via a combination of cylindrical lenses and a distribution optic. The complete system comprises two of these expander units, along with two retro-reflection mirror units and a push beam.

All of the PICAS optical systems are pre-aligned and securely glued in place prior to cold atom testing. This approach ensures that the system is robust, immune to vibration, requires no user-alignment, and can be used in any orientation with respect to gravity.

Product Specifications

Atom Sources	Rubidium (Cs available soon)
Electrical Connections	4 PCB pin connections (2 per source)
Typical Laser Input Requirements	2x 20mW and 0.4mW repump for 2D MOT 1x 2mW and 0.1mW repump for Push Polarization keyed on fast axis of fiber
Optical Input	3x FC/APC Fiber Optic Inputs
Magnetics	Incorporated MAG-2000 magnet assembly
User Adjustment of Optics / Magnets	Via single incorporated XY translation stage
Temperature Range	Vacuum cell up to 225 °C (with PICAS removed)
Dimensions	67mm x 67mm x 93mm

ColdQuanta | 303.440.1284

info@coldquanta.com

To Request a Quote Go to www.coldquanta.com

SOLITON Laser- und Messtechnik GmbH

Talhofstrasse 32 • 82205 Gilching • Telefon +49 (0) 8105-7792-0 • Telefax +49 (0) 8105-7792-77 http://www.soliton-gmbh.de • E-mail: info@soliton-gmbh.de Distributed by:





Photonically Integrated Cold Atom Source - PICAS

Beam-Folding

Prism

Mirror

Product Options

Product numbers Rubidium-rubidium configuration: PICAS-1000-RB Cesium-cesium configuration: PICAS-1000-CS Vigindrical Fiber Connectors

Push Optics

Mechanical Drawing

Distribution Optic

Cold Atom Source Cell



SOLITON Laser- und Messtechnik GmbH Talhofstrasse 32 • 82205 Gilching • Telefon +49 (0) 8105-7792-0 • Telefax +49 (0) 8105-7792-77 http://www.soliton-gmbh.de • E-mail: info@soliton-gmbh.de